

# OSOM®

## Influenza A&B Test

*When rapid, differential diagnosis is essential  
to your patient's health and comfort.*



✓ **ACCURATE, DIFFERENTIAL DIAGNOSIS**

Differentiates between influenza A & B.  
Three distinct lines.

Capable of detecting some animal and avian  
influenza strains.\*

✓ **QC INSIDE**

Includes 2 extra test sticks for external QC testing.

✓ **NASAL SWAB SPECIMENS**


✓ **FAST**

Less than 1 minute hands on time. Results in 10 minutes.

✓ **ROOM TEMPERATURE STORAGE**



**TEST PROCEDURE / RESULTS**



**FOR MORE INFORMATION ON OUR RAPID TEST PRODUCTS, VISIT US AT [www.genzymediagnosics.com](http://www.genzymediagnosics.com)**

\* Although this test has been shown to detect cultured avian influenza viruses, including avian influenza A subtype H5N1 virus, the performance characteristics of this test with specimens from humans infected with H5N1 or other avian influenza viruses are unknown. See package insert for details.

# rapidreader

An E-news Publication of Genzyme Diagnostics

## OSOM® Influenza A&B Test Package Insert Update

Added H1N1 strain (Mexico/4108/2009) to the Analytical Reactivity section

The OSOM Influenza A&B Test has been shown to detect the Mexico/4108/2009 virus in culture isolates. The estimated detectable level for this strain was  $8 \times 10^6$  EID<sub>50</sub>/ml.

Although this test has been shown to detect the 2009 H1N1 virus cultured from a positive human respiratory specimen, the performance characteristics of this device with clinical specimens that are positive for the 2009 H1N1 influenza virus have not been established.

The OSOM Influenza A&B test can distinguish between A and B viruses, but it cannot differentiate influenza subtypes.

For additional information and updates, please refer to the CDC website:

[www.cdc.gov/h1n1flu/guidance/rapid\\_testing.htm](http://www.cdc.gov/h1n1flu/guidance/rapid_testing.htm)

**genzyme**  
Diagnostics

**BOME**

**Sanayi Ürünleri Dış Tic. Ltd.Şti.**

Sokullu Mehmet Paşa Mah. 6. Sokak No : 6/A  
06460 Sokullu/ Çankaya / Ankara / TÜRKİYE  
Tel : (+90 312) 479 33 30 (pbx) Fax: (+90 312) 478 37 18  
www.bome.com.tr E-mail : info@bome.com.tr